

EpiOcular™ Eye Irritation Test (EpiOcular-EIT) for Identification of Materials Not Requiring Classification and Labeling for Eye Irritation or Serious Eye Damage

Yulia Kaluzhny, Helena Kandárová, Laurence d'Argembeau-Thornton, Thoa Truong, Paul Kearney, Patrick Hayden, and Mitchell Klausner



MatTek
CORPORATION

Why We Need Safety Tests for Products




Fig. 1. To depict the need for governmental regulations over cosmetics, the U.S. Food and Drug Administration prepared this poster in 1933 from photographs of a 38-year-old Ohio woman. The patient's first photo was taken an hour before eyelash dyeing, and the other was obtained during the following month showing bilateral staphylococcal corneal ulcers that resulted in vision of light perception. (Reprinted from Lamb RdeF¹⁰⁹ with permission of Farrar and Reinhart.)




Fig. 3. John H. Draize, PhD, the FDA pharmacologist who developed methods for assessing the ocular and dermal toxicity of drugs and cosmetics (Courtesy of the U.S. Food and Drug Administration History Office).

EpiOcular: Reconstructed human Cornea-like Epithelium




A) *in vitro* EpiOcular tissue (OCL-200)




B) rabbit corneal epithelium.

- **Test system:** non-keratinized multi-layered epithelium reconstructed from primary human epidermal keratinocytes.
- **H&E** stained cross-section of the EpiOcular tissue (NHEK) reveals **highly organised basal and supra-basal cell layers** followed by cells which flatten out as the apical surface is approached, similar to *in vivo* corneal tissue.
- **Endpoints usually measured:** cytotoxicity (MTT assay), permeation, release of inflammatory mediators
- **Applicability and limitations:** Applicable to all types of chemicals and also to intensely coloured chemicals (with use of HPLC/UPLC-spectrophotometry)

EpiOcular: Reconstructed human Cornea-like Epithelium



EpiOcular Test Methods



EpiOcular™ Eye Irritation Test (EIT): OECD TG 492



EpiOcular EIT Test Development

- EpiOcular EIT Protocol was developed by MatTek Corporation in 2007, as a reaction on ECVAM's request for a simple and straightforward assay that could be used in the framework of **REACH** and **7th Amendment to the Cosmetic Legislation** (on demand protocol).
- EIT was developed for **classification and labeling substances for regulatory purposes of chemicals including raw cosmetic ingredients**.
- The protocol and its prediction model was designed to discriminate between ocular irritant / corrosive materials (GHS Categories 1 and 2 combined) and those that require no labeling (GHS No Category).
- In total **112 chemicals** were tested during the development of the protocol by MatTek.



COLIPA Pre-Validation of EpiOcular™ Protocol: Technology Transfer - 2007



- Protocol transfer was tested together with IIVS and pre-validation was done in **7 mainly cosmetic oriented laboratories**.
- Independent **in house validation study** with 60 industrial chemicals was performed by **BASF**. Later on, BASF also conducted testing of pesticides formulations using the EIT method and received excellent results.

Detailed information in:

In-house Validation of the EpiOcular™ Eye Irritation Test and its Combination with the Bovine Corneal Opacity and Permeability Test for the Assessment of Ocular Irritation.

Kolle, S.N., Kandárová, H., Wareing, B., van Ravenzwaay, B., Landsiedel, R. (2011). ATLA 39, 365–387, 2011 (60 chemicals).

EpiOcular EIT Test was Adopted by OECD as TG 492 in 2015

OECD/OCDE


492

Adopted:
28 July 2015

OECD GUIDELINE FOR THE TESTING OF CHEMICALS

**Reconstructed human Cornea-like Epithelium (RhCE) test method for identifying chemicals
not requiring classification and labelling for eye irritation or serious eye damage**

EpiOcular™ EIT Protocol




EpiOcular™ EIT Protocol




MTT assay (3hours)

Dead tissue



Viable tissue



Prediction model:

Eye Irritating chemicals are cytotoxic in specific test conditions of the in vitro eye irritation test

If viability > 60%, **Non Irritant**

If viability < 60 %, **Irritant (“no conclusion can be made” within the interpretation of the OECD TG 492)**

OECD TG 492 can be utilized as part of a tiered testing strategy to identify chemicals that do not require classification and labeling.

All compounds predicted as irritants are undergoing further testing (e.g. using BCOP, ICE...) to identify the irritation potency/severity or to exclude false positive prediction.

EpiOcular™ EIT Protocol

jove B N P M BE EN C EN E DE SE

Eye Irritation Test (EIT) for Hazard Identification of Eye Irritating Chemicals using Reconstructed Human Cornea-like Epithelial (RhCE) Tissue Model

Yulia Kaluzhny¹, Helena Kandarová², Laurence d'Argembeau-Thornton¹, Paul Kearney¹, Mitchell Klausner¹
¹MatTek Corporation, ²MatTek In Vitro Life Science Laboratories

Article Downloads Comments Metrics f t e p + 2



0:05 Title
1:43 Reconstructed Human Cornea-like Epithelial (RhCE) Tissue Treatment Preparation (Day 0)
3:08 Pre-treatment and Test Article Exposure (Day 1)
4:42 Rinsing
6:13 MTT Viability Assay (Liquid Test Day 1; Solid Test Day 2)
8:19 Results; Representative EIT Results for 10 Test Articles
9:36 Conclusion

This article is **Open Access**.

Recommend JoVE to Your Librarian

Related Videos

- An In Vitro Skin Irritation Test (SIT)...
Published 7/13/2009
- Laser-induced Chronic Ocular Hypertension Model...
Published 12/04/2007
- Intravitreal Injection for Establishing Ocular...
Published 10/01/2007
- A Laser-induced Mouse Model of Chronic Ocular...
Published 8/14/2013

Translate text to:

Summary

We have developed an eye irritation test which utilizes a three dimensional reconstructed human cornea-like epithelial (RhCE) tissue model. The test is able to discriminate between ocular irritant and corrosive materials (GHS Categories 1 and 2 combined) and those that do not require labeling (GHS No Category).

Date Published: 8/23/2015, **Issue** 102, **doi:** 10.3791/52979

Keywords: Bioengineering, Issue 102, alternative methods, cornea, EpiOcular, EURL ECVAM, eye irritation, in vitro, MTT viability assay, ocular irritation, REACH, reconstructed ocular tissue model, RhCE

<http://www.jove.com/video/52979/eye-irritation-test-eit-for-hazard-identification-eye-irritating>

EpiOcular assays for testing chemicals and formulations – ET-50 protocols




EpiOcular time-to-toxicity (ET-50) protocols

- Protocols exist for more than 25 years
- Suitable for the assessment of the irritation potency, tolerance and mildness of a cosmetic formulations – AVON and Mary Kay use the EpiOcular assay for the mildness testing of their products for more than 20 years!
- Specific protocols were developed for surfactants and surfactant based formulations. Validation was conducted by IIVS and Colgate in 1998. Data from the validation were submitted to ECVAM in 2000 and resubmitted in 2004. Applicability domain increase was requested. New studies initiated within CON4EI project (2015-2017).
- Accepted by EPA for the assessment of Antimicrobial cleaning products. Level of irritation is assessed as tissue viability in MTT assay, determination of ET-50 value.
- Different test design and PM than EpiOcular EIT – more time-points required, but more information can be obtained compared to the EIT.


EpiOcular time-to-toxicity (ET-50) protocols

- **Neat Method:** used for non-water soluble materials or for exposing test articles undiluted.
- **Dilution Method:** applicable to water-soluble materials with a specific gravity of ≥ 0.95 and requires an initial dilution of the test article to 20 % in water. Recommended for surfactants and rinse-off cosmetics.
- **Sub-Draize Mildness Testing:** involves applying neat test articles to the EpiOcular tissue model and is used to differentiate between materials for which standard Rabbit Eye Draize testing is insensitive (eye care cosmetics).

Example of the ET-50 Protocol: Neat Method



Example of the ET-50 Protocol: Neat Method



Draize Score

0-15
15.1-25
25.1-50
50.1-110

Irritancy Classification

Non-irritating, Minimal
Mild
Moderate
Severe, Extreme

Example

PEG-75 Lanolin, Tween 20
Pareth 25-12
1% Triton X-100
5% Benzalkonium Chloride

EpiOcular ET-50 (min)




>60
30-60
3-29.99
<3

Example of the ET-50 Protocol: Neat Method

Assessment of Ocular Irritation Ranges of Market-Leading Cosmetic and Personal-Care Products Using an *In Vitro* Tissue Equivalent – SOT 2002 – The Toxicologist 66, 243, 2002. N E McCain, R R Binetti, S D Gettings, and B C Jones, **Avon Products**, Inc., Suffern, NY, USA

Formulation Category	Number of Samples	Mean ET50 (min)	Median ET50 (min)	Range of ET50 Scores (min)
Surfactant Formulas				
Adult Shampoo	9	27.0	28.8	17.1-35.1
Children's Shampoo	7	58.2	45.4	37.9-135.4
Baby Shampoo	4	74.9	80.3	31.9-107.1
Baby Wash	4	114.5	112.6	76.9-155.9
Eye Area Cosmetics				
Liquid Eyeliner	15	532.2	127.0	64-1440
Eye Cream	11	581.8	454.0	102-1200
Mascara	20	655.1	626.5	199-1080
Pencil Eyeliner	3	1440.0	1440.0	1440
Creams/Lotions				
Adult Lotion/Cream	17	511.4	440.8	234-1440
Sunscreen	15	673.6	486.0	252-1440
Baby Lotion	2	1200.0	not applicable	1200
Non-Surfactant Haircare				
Children's Styling Gel	2	172.7	not applicable	126-219.3
Adult Conditioner	5	279.6	256.0	100-530
Children's Conditioner	2	728.0	not applicable	256-1200
Children's Hair Detangler	3	868.4	1200.0	205.1-1200

Example of the ET-50 Protocol: Neat Method



All of the formulations tested in this study, with the exception of six adult shampoos, had individual ET₅₀ scores greater than 31 minutes and are classified as mild to minimal/non-irritating. The six adult shampoos with ET₅₀ values less than 31 minutes are classified as moderate irritants.

Summary

- **EpiOcular EIT test protocol has been adopted on July 28, 2015** as the new **OECD TG 492**. This test enables reliable identification of substances that do not require labelling.
 - If a chemical classifies as NI in the EIT, no further testing is required.
 - In case of positive results from the EIT it is necessary to conduct further *in vitro* testing to exclude corrosive nature of the test article or false positive prediction. Such a test could be BCOP, ICE or STE if the applicability domain is not limited.
- AVON, Mary Kay and some other **major cosmetic companies have established a reference database** of *in vitro* ocular irritation using the **EpiOcular ET-50 method** for a cross-section of marketed cosmetic and personal care products.
- The *in vitro* EpiOcular tissue model used in the ET-50 test design is able to distinguish differences between moderate **“mild”** and **“ultra-mild”** formulations.

Training in performing the OECD TG 491 and the OECD TG 492 in Brazil

Project Leader: Dr. Talita Miguel Marin



The Brazilian National Alternative
Methods Network (RENAMA)



Ministério da
Ciência, Tecnologia,
Inovações e Comunicações



PLATAFORMA REGIONAL DE MÉTODOS ALTERNATIVOS AO USO DE ANIMAIS DE EXPERIMENTAÇÃO DO MERCOSUL



CURSO

Inscrições até 01 de outubro de 2018

**Métodos de avaliação do potencial de irritação e
corrosão ocular que utilizam modelos *in vitro* de
córnea (OECD TG 491 e 492),**

Coordenadora: Dra. Talita Miguel Marin

Data: 26/11/18 à 30/11/18.

Local: Laboratório Nacional de Biociências – LNBio/ CNPEM
Rua Giuseppe Máximo Scoffaro 10.000, Cidade Universitária,
Campinas, São Paulo.

Inscrições

(Brasil):

<https://goo.gl/forms/JOAhtv4BapFBn8rm1>

Conteúdo: capacitar os participantes para a realização do ensaio TG OECD 491 para identificar substâncias que causam danos severos oculares e substâncias que não produzem dano. O curso também abordará o ensaio baseado em córnea humana reconstituída (TG OECD 492) e as novas metodologias envolvidas para emular a resposta humana (human-on-a-chip).

Público-alvo: Estudantes de pós-graduação, profissionais e pesquisadores que atuem na área de Métodos Alternativos e/ou experimentação animal nos países do MERCOSUL tanto do setor público quanto privado.

BRAZILIAN CENTER FOR RESEARCH IN ENERGY AND MATERIALS – CNPEM Campinas, São Paulo



Project Leader: Dr. Talita Miguel Marin

Thank you for your attention!

For more information please contact:



Dr. Yulia Kaluzhny:
ykaluzhny@mattek.com



Dr. Helena Kandarova:
hkandarova@mattek.com

Dr. Patrick Hayden:
phayden@mattek.com



MatTek | IN VITRO
LIFE SCIENCE
LABORATORIES